


AUTO STERILIZER AID

BioEra's Auto Sterilizer Aid utilizes infrared heat produced by ceramic core element for complete sterilization. The ceramic core element ensures absolute sterilization within 5-7 seconds at a temperature of 1500°F (815°C).



SAY NO TO Flame, Burners, Gas, Spirit Lamps
SAY YES TO Sterilize with Handy, Small Foot-Print, Working on Power-Input, Speedy & Secured Sterilization with BioEra's Auto Sterilizer Aid

CAT. NO.: BE/ASA/01



FEATURES

- Hassle free complete sterilization within 5-7 seconds at a temperature of 1500°F (815°C)
- Utilization of ceramic core element to ensure absolute sterilization
- Sterilization based on infrared heat
- Adjustable angle of heater to ensure user convenience for sterilizing various surface area with different heights (like tubes, flasks, wire loops etc.)
- Heater angle range adjustable up to 55 degree angle
- Length of sterilizer heater : 150mm
- An ideal sterilization system for both aerobic and anaerobic work areas
- Power requirement : AC220V; 50/60Hz
- Net weight : 1.2kg
- Dimensions (D X W X H) : 159 X 100 X 260 mm