

# STEAM STERILIZERS (AUTOCLAVES) : MODEL STANCE

BioEra's Steam Sterilizers are made to suffice the European standard EN 13060. This system offers a perfect solution for sterilization and multiple utility in various application areas.



MODEL: STANCE [18 Liters]  
CAT. NO.: BE/CI/BSS/SA-18

## FEATURES

- Class N autoclave (according to European Standard EN 13060)
- LED Display
- High precision temperature control of 0.01°C
- Automatic steam venting
- Overheat protection system will cut-off power supply in case of overheating
- Pressure monitoring through pressure gauge
- Door protection locking system prevents the door from opening till presence of pressure in the chamber
- Water drainage through water drainage valve
- Power Requirements – AC 220 V – 240 V; 50 - 60Hz
- Overall Dimensions – 445 X 400 X 605 mm
- Net Weight – 30kg

## TECHNICAL SPECIFICATION

FEATURES	DESCRIPTION
Model	STANCE
Cat. No.	BE/CI/BSS/SA-18
Chamber Volume	18 Liters
Sterilization Temperature	121°C & 134°C
Display	LED
Weight	30kg