

DISTILLATION UNIT MODEL FRILL

BioEra's distillation unit model frill is designed as a miniature system to serve laboratory distill water requirements. The machine comes with compact design and easy to operate technology. The unit has capacity to generate 4 liters of distill water with 1.5 liters per hour output. Machine comes with fan cooling system with which steam gets condensed and is collected in water collection chamber. As the water is not wasted in cooling the steam it is eco-friendly. The machine has auto-shutdown mechanism when the process is completed which safeguards the machine.



No Water Wastage



4 Ltr Capacity



In Built Cooling Fan



Water Distillation Capacity
1.4 Ltr per Hour

CAT. NO.: PDU-310360

FEATURES

- LED indicator
- Energy efficient & Eco-friendly
- Machine auto shuts after the water is over in the chamber
- Comes with collection chamber for easy storage & portability
- The distillation unit is essential laboratory aid with easy to operate & maintain process
- Exterior of machine is made up of engineering plastic & interior chamber is made up of stainless steel 304
- Water tank capacity 4 liter
- Water distillation capacity- 1.5 Liter/Hour
- No wastage of water because of continuous vaporization & condensation during distillation cycle
- In-Built fan cooling system for condensation
- No Need of tap water connection
- Comes with detergent for washing of residue



Collection Chamber



Fan Cooling System