

EONPIPETTE (AUTOCLAVABLE)

FEATURES SET - COMMON

- Best quality raw material
- Extremely durable
- Designed for exceptional user ease
- Ensures high precision, accurate and consistent results
- Soft touch tip ejection
- Easy adjustment within its volume range
- Tipcones are designed to fit most pipette tips
- Non-slippery firm grip
- A handy volume setting dial for single handed operation
- Easy to calibrate and maintain
- Tipcone and Tipsejector are resistant to most of the Laboratory Chemicals
- Spare parts and Accessories are readily available
- One year conditional warranty
- EonPipette Model Grand & Model Eco are Fully Autoclavable
- One set of Variable Volume EonPipette comprises five pipettes covering the widest volume range from 0.1 μ l to 5000 μ l

EONPIPETTE –MODEL GRAND
(FULLY AUTOCLAVABLE)



FULLY
AUTOCLAVABLE
Universal Tipcone to fit most
pipette tips
Covers VOLUME RANGE from
0.1 to 5000 μ l



EONPIPETTE – MODEL ECO
(FULLY AUTOCLAVABLE)

WARRANTY

BioEra's EonPipettes has been manufactured and tested under stringent quality control to ensure accuracy and precision. These units are warranted to be free from defect in material and manufacturing process for a period of one year from the date of purchase. The warranty against defect in manufacturing process or material will be accepted as long as the product was used in compliance with the instruction for maintenance and use. Pipettes exposed to any physical or chemical damage will not be accepted under any warranty terms. Calibration and servicing charges are not included under warranty terms.

VARIABLE VOLUME EONPIPETTES MODEL ECO

CAT. NO.	DESCRIPTION	ACCURACY (%)	PRECISION (CV %)
VV0005	EonPipette Variable Volume ECO Range 0.1 - 2 μ l	\pm 12.0-1.5	\leq 6.00 – 0.70
VV0010	EonPipette Variable Volume ECO Range 2 - 20 μ l	\pm 2.5 – 1.0	\leq 1.50 – 0.30
VV0020	EonPipette Variable Volume ECO Range 20 - 200 μ l	\pm 1.8 – 0.6	\leq 0.50 -0.15
VV0100	EonPipette Variable Volume ECO Range 100 - 1000 μ l	\pm 1.5 – 0.6	\leq 0.30 – 0.15
VV0200	EonPipette Variable Volume ECO Range 500 – 5000 μ l	\pm 1.2 – 0.6	\leq 0.30 – 0.15

VARIABLE VOLUME EONPIPETTES MODEL GRAND

CAT. NO.	DESCRIPTION	ACCURACY (%)	PRECISION (CV %)
VV23001	EonPipette Variable Volume GRAND Range 0.1 - 2 μ l	\pm 12.0-1.5	\leq 6.00 – 0.70
VV23002	EonPipette Variable Volume GRAND Range 2 - 20 μ l	\pm 2.5 – 1.0	\leq 1.50 – 0.30
VV23020	EonPipette Variable Volume GRAND Range 20 - 200 μ l	\pm 1.8 – 0.6	\leq 0.50 -0.15
VV23100	EonPipette Variable Volume GRAND Range 100 - 1000 μ l	\pm 1.5 – 0.6	\leq 0.30 – 0.15
VV23500	EonPipette Variable Volume GRAND Range 500 – 5000 μ l	\pm 1.2 – 0.6	\leq 0.30 – 0.15

SCHEME Purchase

Ten Variable Volume EonPipette Model Grand

Receive Free

Four Variable Volume EonPipette Model Grand of your choice

VERTICAL FLAT PIPETTE RACK

FEATURES

- Capacity to occupy Five Pipettes
- Made up of high quality acrylic material
- Durable
- Ergonomic design
- Attractive look
- Compact Size
- Perfect to be used with most of the Brand of Pipettes available in market
- One year conditional warranty



CAT. NO.	DESCRIPTION
VR 001	Vertical Flat Pipette Rack