

MAGNETIC STIRRER MODEL EARLE

(Available with Heating Function)

BioEra's Magnetic Stirrer Model Earle is based on Electromagnetism Driven Mechanism and Nano Heating Technology for unmatched performance. The system comes with heating function.



Comes with
free
Magnetic Needle,
Support Clamp &
Temperature Sensor

CAT. NO.: BE/MST/H/03

Comes with Temperature Sensor & Support Clamp

FEATURES

- Unmatched aesthetics and advanced design
- Electromagnetism driven mechanism
- Maintenance free system
- No risk of demagnetization
- User friendly knobs
- LED display provides easy reading of the speed and time
- Ceramic working platform is easy to clean, waterproof and resistant to strong acids & alkali
- Automatic reverse rotation every 30 seconds provides excellent mixing



TECHNICAL SPECIFICATION

FEATURES	SPECIFICATIONS
Model	Earle
Cat. No.	BE/MST/H/03
Stirring Positions	1
Working Platform Material	Ceramics
Speed Range	15 - 1500
Display	LED
Maximum Stirring Quantity	5000ml
Stirring Mechanism	Electromagnetism Driven
Working Platform Size	5 inch
Reverse Rotation	Automatic after every 30 Seconds
Temperature Range	50 – 310°C
Temperature Accuracy	±0.1°C
Power Requirement	AC220V; 50/60Hz
Dimensions (W x D x H) mm	232 x 45 x 312
Net weight	1.7kg