


MAGNETIC STIRRER MODEL ADROIT

BioEra's Magnetic Stirrer Model Adroit is made with advanced magnetic coil technology to deliver efficiently. Touch key function and ultrathin design makes this model ideal to enhance aesthetics of laboratory. In addition the resistance to corrosion is offered due to platform made up of tempered glass.



One of the most Ultrathin design in World
 LED Display with Touch Key Operation
 Glass Ceramic Platform
 Auto Reverse Rotation Function
 Maximum Stirring Quantity – 5000ml

CAT. NO.: BE/MST/02

FEATURES

- Unmatched aesthetics and advanced design
- Electromagnetism Driven Mechanism
- Maintenance free system
- No risk of demagnetization
- Touch key function
- LED display provides easy reading of the speed and time
- Glass ceramic working platform is easy to clean, waterproof and resistant to strong acids & alkali
- Automatic reverse rotation every 30 seconds provides excellent mixing



TECHNICAL SPECIFICATION

FEATURES	SPECIFICATIONS
	MODEL ADROIT WITHOUT HEATING FUNCTION
Cat. No.	BE/MST/02
Stirring Positions	1
Working Platform Material	Glass Ceramics
Speed Range	15 - 1500
Display	LED
Maximum Stirring Quantity	5000ml
Stirring Mechanism	Electromagnetism Driven
Working Platform Size	5 inch
Reverse Rotation	Automatic after every 30 Seconds
Power Requirement	AC220V; 50/60Hz
Dimensions (W x D x H)	232 x 45 x 312 mm
Net weight	1.7kg