

# PROGRAMMABLE MICROCENTRIFUGE: MODEL FIST

BioEra's Programmable Microcentrifuge is powered with attractive design and sophisticated technology. This model is designed to present durability, reliability and safety to the user. Features like timer, safety interlock switch, monitor window etc. boost up this unit as a perfect system to employ for all volume spinning needs. In addition this model is boosted with the features of attractive look and stunning color.



CAT NO.	DESCRIPTION
PMC-260606	Programmable Microcentrifuge Model Fist (Angle Rotor) (300 – 4000 rpm)
PMC-260607	Programmable Microcentrifuge Model Fist (Microplate Rotor) (500 – 2000 rpm)

Complementary 50ml Tube Rotor & Adaptors For 0.5ml / 1.5/2ml / 5ml / 7ml / 15ml  
(Adaptors 2 No. Each)

## FEATURES

- Unique design and attractive look
- Super low working noise i.e.  $\leq 48\text{dB}$
- Easy to handle with its sleek design and low weight
- Fiberglass plastic body ensures durability reliability and low weight
- High intension PC lid ensures safety at high speed working of instrument
- Safety interlock switch- actuated by unit lid
- Tube position on the rotor are numbered for easy identification
- Aid of rubber grommets at the base of centrifuge unit for the better gripping to surface
- Brushless Rotor
- Timer: 15 seconds to 99 minutes
- Additional accessories like rotor and adaptors available at an additional cost
- Rotor capacity: 8 places X 0.5ml / 1.5ml / 2ml / 5ml / 7ml / 10ml / 15ml / 50 ml
- Centrifugation speed: 300 – 4000 rpm
- Maximum RCF: 2460 g
- Additional accessories like adaptors available for 0.5ml / 1.5ml / 5ml / 7ml / 15ml tubes
- Power Supply: 110 V – 240 V; 50 Hz
- Dimensions: 370 X 330 X 175 mm
- Approximate weight: 11 kg

## ACCESSORIES

CAT NO.	DESCRIPTION	PACK SIZE
MC-ACC-09	Set of Tube Adaptors 1.5ml / 5ml / 7ml / 15ml	2 Nos Each