


ROTARY VANE VACUUM PUMP

Rotary vane vacuum pump is connected with shaft and provides high speed. It is smaller in size, compact and convenient to use.



CAT. NO.: BE/CI/RE-03

FEATURES

- Table Top Design, compact in size saving space and easy to move
- Elegant design and portable to handle
- Pumping rate : 0.25~30L/S
- Gas type : Clean and dry air without other mixture at room temperature
- Single phase asynchronous motor
- Easy Oil drain port
- Large window to observe the oil level at any time
- Water-proof switch, safe and reliable
- Air bleed speed m^3/h (L/S) : 21.6
- Extreme pressure (Pa) : $\leq 4 \times 10^{-2}$
- Rotating speed (r/min) : 1400
- Motor power (kw) : 0.75
- Noise (dBA) : 68
- Oil volume (L) : 2
- Size mm : 565 x 200 x 342
- Weight (Kg) : 15kg

