


VORTEX MIXER

BioEra's Vortex Mixer is next generation model with smallest foot-print. It is handy and most efficient device for mixing. The instrument is light in weight with specially designed foot pad to hold the Vortex Mixer firmly in place when functioning.



CAT. NO.: BE/VM/01

Touch function



FEATURES

- Small foot-print and compact in size
- Touch function
- Orbital movement with Fixed speed of 3000rpm
- Orbital Diameter of 5.5mm
- Specially designed foot pad can hold the vortex mixer firmly in place when in mixing
- Complete mixing within 5 seconds for full volume of 50ml
- Power requirement : AC220V; 50/60Hz
- Net weight : 998gm
- Dimensions (Φ x H) : 100 X 78 mm