

# pH METER MODEL ADEL

BioEra's pH meter is unique economical and attractive. Body is made up of engineering plastic with Large LCD. The machine is equipped with temperature sensor for precise and accurate measurement of pH during temperature fluctuations.



Made in Plastic Body  
with Attractive Design  
Big LCD with Backlight  
Measures pH, mV  
& Temperature

CAT. NO.: BE/PHTT/03

## FEATURES

- Attractive design with body made of plastic
- Three parameters measurement pH, mV & Temperature
- Big LCD with Illuminous back light display
- pH Range : 0.00~14.00 ; pH Accuracy :  $\pm 0.01 \pm 1$  bit ; pH Resolution : 0.01
- mV Range :  $\pm 1999$  ; mV Accuracy :  $\pm 0.1$  ; mV Resolution : 1
- Adjustable probe holding stand
- Automatic temperature compensation
- Dimensions (W x D x H) - 160 x 190 x 70 mm

## TECHNICAL SPECIFICATION

FEATURES	DESCRIPTION
Cat. No.	BE/PHTT/03
Measurement Parameters	pH / mV (ORP)/ °C
Range	pH : 0.00 ~ 14.00 mV : $\pm 1999$
Resolution	pH : 0.01 mV : 1
Accuracy	pH : $\pm 0.01 \pm 1$ bit mV : $\pm 0.1\% \pm 1$ bit
Calibration	Auto Calibration with Slope Display
Power Requirement	Dc 9v Power Adapter
Dimensions (W x D x H) mm	160 x 190 x 70
Net weight	0.75 Kg